COMBILITE® , COMBISTAR® SERIES

GCE COMBINATION VALVES FOR USE WITH MEDICAL GASES

EN
DE
FR
NL

INSTRUCTIONS FOR USE

Document Nr: 735800000063
Date of issue: 01/03/2010
Revision Nr: 8

COMBILITE®, COMBISTAR® SERIES

GCE COMBINATION VALVES FOR USE WITH MEDICAL GASES

1. Foreword

The product complies with essential requirements of 93/42/EEC Medical Device Directive and 99/36/EEC Transportable Pressure Equipment Directive.
The valve is designed according to EN ISO 10524-3 standard.

2. Intended use

Combilite, Combistar Combination Valves are designed to be fitted to gas cylinders used for medical gases. These combination valves with a gas cylinder form medical gas packages used either as gas supply point for medical devices (anaesthetic devices, ventilating devices, incubators etc.) or for direct gas supplying to a patient's breathing mask or cannula.

GCE Combination Valves are intended to be used with the following medical gases:

- | | |
|---|--|
| <ul style="list-style-type: none"> - oxygen; - nitrous oxide; - air for breathing; - helium; - carbon dioxide; | <ul style="list-style-type: none"> - xenon; - specified mixtures of the gases listed; - air for driving surgical tools; - nitrogen for driving surgical tools. |
|---|--|

3. Operational, transport and storage safety requirements

	Operations	Transport	Storage
<div style="display: flex; align-items: flex-start;"> <div style="margin-right: 10px;"></div> <ul style="list-style-type: none"> • Keep away the product and its associated equipment from <ul style="list-style-type: none"> - heat sources (fire, cigarettes, ...), - flammable materials, - oil or grease (take a great care in the use of hand creams) - water, - dust. </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ ✓ ✓ ✓ ✓ </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ ✓ ✓ ✓ ✓ </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ ✓ ✓ ✓ ✓ </div>
<div style="display: flex; align-items: flex-start;"> <div style="margin-right: 10px;"></div> <ul style="list-style-type: none"> • The product and its associated equipment must be prevented from tipping over, turning over or falling. </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ </div>	<div style="display: flex; flex-direction: column; gap: 5px;"> ✓ </div>
<div style="display: flex; align-items: flex-start;"> <div style="margin-right: 10px;"></div> <ul style="list-style-type: none"> • Maintain always oxygen cleanliness standards. </div>			

- Use only the product and its associated equipments in well ventilated area.

Before the first use the product shall be located in its genuine package. In case of removal from service (for transport, storage) GCE recommends using the genuine package (including inner cushioning materials).

National laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

Operating conditions

 -20/+60 °C

 10/100 %

 600/1200 mbar

Storage and transport conditions

 -30/+60 °C

 10/100 %

 600/1200 mbar

- In case of valve storage at temperature below -20°C do not use the valve until its temperature reaches at least -20°C.

- For the valves designed to be used with mixture of gases O₂+N₂O, the lowest operating temperature is +5°C. In common use of the valve, frosting can appear on the valve surface, which is caused by common physical reaction inside the valve when high pressure in the valve is reduced to low pressure (Joule-Thomson effect). Check that all associated equipment is connected to the valve via a hose of the length of at least 2 metres.

4. Personnel instructions and training

According to the Medical Device Directive 93/42/EEC, the product owner must ensure that all users and persons handling the product are provided with operating instructions & performance data of the product.

The product owner shall ensure that all users and persons handling the product are fully trained to carry out that operation. Trainees need to be supervised by an experienced person.


Contact GCE if further information is required on GCE product training program.

5. Product description

Combination valve combines the function of shut-off valve of a high-gas cylinder and pressure regulator for the use with medical gases. Gas from the cylinder is first controlled by the main shut-off valve and then passes through the pressure regulator and finally delivered to the patient through the flow outlet or the pressure outlet. Outlet positive pressure is fixed by the manufacturer and each valve is provided with (external or internal) low-pressure relief valve against undesirable exceeding of outlet positive pressure.

There are three basic alternatives:

- valve outlet via quick-coupler, the outlet flow is a constant one.
- valve outlet via calibrated nozzles, outlet flows can be changed by the control knob
- the combination of both alternatives

(Typical combination valve configuration)

A – Inlet stem

The valve is fitted to gas cylinder by a threaded inlet stem. The inlet stem can be taper threaded or parallel threaded of different size depending upon cylinder size and material.

B - Filling port

Filling port is provided for filling the gas cylinder. It includes a non-return valve (NRV). The NRV means that special filling adaptors are required to vent gas from the cylinder during the filling process.

C – Pressure gauge

A cylinder contents pressure gauge is provided. This can be active or non-active. Versions fitted with active pressure gauge enable to check cylinder content at shut-off valve in both ON and OFF positions. The versions fitted with non-active pressure gauge enable to check cylinder content at shut-off valve in ON position only.

D – Shut-Off valve

The combination valve is provided with a shut-off valve to isolate the gas from the cylinder from the other valve functions. It must be turned on during cylinder filling and patient therapy.

E – Residual Pressure Valve

GCE combination valve contains a residual pressure valve whose function is to retain a minimum positive pressure in the gas cylinder to avoid contamination of the cylinder content by atmospheric air.

F, G – Flow control head “F” and Flow outlet “G” (option)

GCE combination valves can be delivered with a flow control head “F”. This function is used to supply gas flow rates (l/min) at atmospheric pressure directly to a patient through the flow outlet “G”, for instance through a cannula or a facemask.

The flow outlet “G” is fitted with hose fitting (hose nipple) or a threaded type (for humidifier for instance) or with combination of both types..

H – Pressure outlet or quick-coupler (option)

The combination valve may be fitted with a pressure outlet. The pressure outlet is the outlet direct from the low-pressure part of the valve and it is fitted with a gas specific medical quick connector also called “quick coupler” (See appendix Nr 2). The user can connect another piece of equipment to this outlet with a gas specific male probe. The quick connector self seals when the male probe is disconnected. This outlet is for supplying gas at a controlled pressure to power medical devices, for instance medical ventilator.

Excess flow device (option)

Excess flow safety valve can be provided and is located in the valve inlet stem. This provides extra safety in case of the valve breaking above the valve inlet stem (for instance at cylinder falling). It has no effect on valve function during normal service

Note: Product colour (especially guard, flow control knob and shut-off valve) might not follow any gas colour coding.

6a.1 Guarding

Valve guarding should be provided where practicable. For cylinders of 5 litres and above a guard must be provided. It shall meet the drop test requirement according to ISO11117 carried out by the gas cylinder owner. It should be designed to enable an easy access to all controllers of combination valves.

- **Installer must ensure that the total medical gas package weight filled with its medical gas does not exceed the maximum weight specified for respective guard (the weight of drop test tested sample). The maximum weight specified for the guard may vary if the guard is fitted or not with a bed hanger.**

- **Installer shall ensure that the guard does not prevent correct valve function, specifically the relief valve or hinder the easy reading of labels.**

6a.2 Fitting to the gas cylinder

Visual Inspection before fitting to the gas cylinder

- Check if there is visible damage on combination valve and gas cylinder (especially inlet stem thread of the valve and thread of gas cylinder). If it shows signs of external damage, remove the valve or gas cylinder from service and suitably identify its status.
- Visually check if combination valve or gas cylinder is contaminated; if needed, apply the cleaning procedure hereunder (for gas cylinder, refer to the cylinder manufacturer cleaning procedure).
- Check if the service due or the total life time of GCE product and the gas cylinder has not been exceeded – see Chapter 9. If service or total life time has been exceeded, remove GCE product from service & suitably identify its status.
- Ensure that the product inlet stem is compatible with the medical gas cylinder (gas/ thread type).
- Check the presence & the integrity of seals / sealing materials (PTFE tape, O-ring) on inlet stem.
- Gas cylinder pressure rating must be at least equal to or greater than the rated combination valve inlet pressure.

- **The product is designed only for use with the gas specified in its labelling. Never try to fit it on gas cylinder for a different gas.**

Fitting of valve to gas cylinder

- Observe oxygen cleanliness standards during fitting.
- In fitting of combination valve to the cylinder follow torque figures specified in EN ISO 13341.

- **Do not exceed the maximum torque allowed for the cylinder. Too high torque may damage the valve thread preventing reuse. Certain types of cylinders may require a torque lower than stipulated in EN ISO 13341. Check torque value with your cylinder supplier.**

- **Tighten inlet stem thread only contacting square surfaces on combination valve body.**

- **Multiple installation or removal can damage inlet stem thread.**

- **Take care not to damage external components for instance pressure gauge.**

6a.3 Filling operations

Visual inspection before filling

- **Remove protective caps from filling port, pressure and flow outlet (if any). Keep the caps in a safety place for reuse in transport or storage.**

- Check if there is visible damage on the product, its parts and filling adaptor. If it shows signs of external damages, remove it from service and suitably identify its status.
- Visually check if the product, its parts, especially filling port and filling adaptor is contaminated; if needed, apply the cleaning procedure hereunder (For the medical gas cylinder and the filling adaptor, refer to the cylinder manufacturer for cleaning procedure if needed).

- Check if the service or the total life time of GCE product and gas cylinder has not been exceeded – see Chapter 9. If service or life time has been exceeded, remove the product from service & suitably identify its status.
- Ensure that the combination valve filling port is compatible with the filling adaptor (gas/type).
- The filling pressure (settled pressure at 15°C) shall not exceed rated combination valve inlet pressure.

Filling adaptor connection

- Observe oxygen cleanliness standards during filling.

- **The combination valve must be vented prior to connecting the filling adaptor.**

- Combination valve internal venting

- **Venting through the combination valve quick-coupler must only be performed with appropriate quick-coupler probe. Contact GCE for further information about quick-coupler probe available for your product.**

For combination valves fitted with flow control head:

- ensure that the shut-off valve is turned off,
- set flow control knob on the maximum position - ensure that the flow control knob engages correctly,
- wait until the combination valve has completely been vented,
- set flow control knob on ZERO position - ensure the flow control knob engages correctly,

For combination valves without flow head:

- ensure that shut-off valve is turned off,
- insert a quick-coupler probe in the combination valve quick-coupler (Refer to appendix Nr 2 to get information about connecting / disconnecting procedure of quick coupler type on your product).
- wait until the combination valve has completely been vented,
- disconnect quick-coupler probe,

- Only use genuine GCE filling adaptors. Contact GCE to get information about available filling adaptor for your product.
- Connect genuine GCE filling adaptor. Do not use excessive torque when connecting the filling adaptor.

- **GCE do not recommend using power drivers or impact drivers to connect filling adapter (excessive force & shocks could result in damage to the filling port). Contact GCE to ensure that the filling port features on your product are compatible with power drivers or impact drivers.**

- **Non GCE filling adaptors can damage the non-return valve in the filling port and impact the safety of the product.**

Medical Gas Package Filling

If the combination valve is used to dispense medical gases within the EC, only gases which have been manufactured according to the requirement of Good Manufacturing Practice (GMP), detailed in the EC Good Manufacturing Practice guide and defined in EC Directive 91/356EEC Annex 6 should be used. Further guidance can be found in EIGA document IGC Doc 99/03/E Good Manufacturing Practice Guide for Medical Gases.

- Prior to filling carry out checks (refer to appendix nr. 3 “Pre / Post filling checks”).
- During filling, the shut-off valve must be fully open.
- After filling ensure that shut-off valve is closed. When turning off at the shut-off valve do not use excessive force (Maximum recommended closing torque is 3 Nm).
- The filling line (including the filling high-pressure hose) must be vented to disconnect the filling adaptor).
- Before to release the medical gas primary package for patient use, perform functional and leak-tight checks Refer to appendix nr 3 “Pre/Post filling checks”, which includes leak-tight test according EIGA document IGC Doc 99/03/E Good Manufacturing Practice Guide.

A special attention should be paid to the following items:

- pressure gauge functionality,
- for combination valves fitted with a flow head, check that there is a gas flow at each setting (for instance, by sound or special leak detector),
- integrity of product label, gas label, anti-tamper label (depending on the type), flow marks on flow control head,
- presence of protective caps & plugs.

- **Do not use leak detection fluids containing ammonia or other substances which can promote stress corrosion cracking in brass materials or can affect other materials. On oxygen equipment use only leak detection fluids approved for use with oxygen systems. Contact GCE if further information is needed on approved leak fluids for your product.**

- **If the relief valve opens, return the valve to GCE for service.**

- After completion of the functional and leak-tight checks, refit the filling port, pressure outlet and flow outlet protection cap as soon as possible. Before refitting the caps ensure that they are clean. Use only genuine GCE caps. If required, check that the anti-tamper label is in position.

6a.4 Removal from cylinder

- Observe oxygen cleanliness standards during removal.
- Cylinder venting – follow Chapter 6a.5
- Combination valve internal venting
 - For combination valves fitted with flow control head:
 - ensure that the shut-off valve is turned off,
 - set flow control knob on the maximum position - ensure that the flow control knob engages correctly,
 - wait until the combination valve has completely been vented,
 - set flow control knob on ZERO position - ensure that the flow control knob engages correctly,
 - For combination valves without flow head:
 - ensure the shut-off valve is turned off,
 - insert a quick-coupler probe in the combination valve quick-coupler (Refer to appendix Nr 2 to get information about connecting / disconnecting procedure of quick coupler type on your product),
 - wait until the combination valve has completely been vented,
 - remove the quick-coupler male probe,
- Connect a GCE genuine filling adaptor. Do not use excessive torque when connecting the filling adaptor.

- **For venting of residual pressure through the filling port use only GCE genuine filling adaptor. Contact GCE for further information about genuine filling adaptor available for your product.**

- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Wait until the whole medical gas primary package has completely been vented.
- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).

6a.5 Cylinder venting

Combination valve fitted with flow control head:

- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Set flow control knob on the maximum position - ensure that flow control knob engages correctly,
- Wait until the cylinder has completely been vented – at least until pressure indicator enters red zone.

- **Perform cylinder venting in well-ventilated areas to avoid oxygen saturation of ambient atmosphere, which could result in an increased danger of fire.**

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).

Combination valve without flow control head:

- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Insert a quick-coupler probe in the combination valve quick-coupler (Refer to appendix Nr 2 to get information about connecting / disconnecting procedure of quick coupler type on your product),

- **Venting through the combination valve quick-coupler must only be performed with appropriate quick-coupler probe. Contact GCE for further information about quick-coupler probe available for your product.**

- Wait until the cylinder has completely been vented – at least until pressure indicator enters red zone

- **Perform cylinder venting in well-ventilated areas to avoid oxygen saturation of ambient atmosphere, which could result in an increased danger of fire.**

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).
- Remove the quick-coupler probe,

- **If in filling process the content of medical gas primary package is quickly vented, GCE recommends that the owner should check that venting parameters (pressure, flow, time, ...) are compatible with the product. Contact GCE for further information about correct venting parameter of your product.**

6b. Distributor operations

6b.1 Before use

- **The hereunder operations are performed for instance by the owner (operator) of the device, trained hospital staff or emergency care staff or trained supplier ensuring delivery of complete medical package to the user, patient.**

Visual inspection before use

- Check if there is visible damage on the medical gas primary package (including product labels and marking). If it shows signs of external damages, remove the product from service and suitably identify its status.
- Visually check if the medical gas primary package is contaminated; if needed, apply the cleaning procedure hereunder.
- Ensure turning-off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).
- Set flow control knob on ZERO position - ensure that the flow control knob engages correctly.

- **Remove caps from flow outlet and/or pressure outlet. Keep caps in a safe place for reuse during transport or storage.**

Leak-tight check before use

- Set flow control knob (if any) on ZERO position - ensure the flow control knob engages correctly.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Visually and by sound check possible leakages from:
 - the filling port (if possible),
 - the flow outlet (if any),
 - the pressure outlet (if any),
 - the relief valve venting hole(s).

- **If any leakage is detected, apply the procedure described in Chapter 6b.3 and return the valve to GCE for service.**

Functional check before use

- Set flow control knob (if any) on ZERO position - ensure that flow control knob engages correctly.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Check if pressure is indicated on pressure gauge. If pressure indicator enters red zone, return the medical package (cylinder with valve) for its refilling.
- Check that there is a gas flow at each set positions in both clockwise and anticlockwise turning direction (for instance by sound or checking presence of bubbles in a humidifier).
- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).
- Reset flow control knob on ZERO position when it's not venting anymore - until the flow control knob engages.
- For combination valve fitted with pressure outlet, ensure it is in working condition by connecting and disconnecting quick-coupler probe (refer to the section hereunder for pressure outlet connection).

6b.2 User outlet(s) connection

List of recognized accessories

To be connected to the flow outlet:

- humidifier extension (the outlet hose nipple is in horizontal position, a "L-shaped" extension may be necessary to use a humidifier, ask for special extension & humidifiers from GCE),
- oxygen tubing connected with breathing mask or cannula,
- gas saver,
- atomizer.

To be connected to the pressure outlet:

- low-pressure hoses,
- flow meters,
- Venturi suction ejectors.

- **Before connecting any accessory or medical device to the regulator, always check that they are fully compatible with the product connection features & the product performances.**

Pressure outlet

- Ensure that the quick-coupler probe is compatible with the pressure outlet features.
- Connect the quick-coupler probe.
- Ensure that quick-coupler probe is well engaged.

- **If the pressure indicator enters red zone and the gas is supplied via pressure outlet, outlet pressure P2 can be lower than 3.6 bar (see appendix nr. 1 for each product).**

- **If pressure outlet is to be connected to medical device that requires high gas flow (for instance pulmonary ventilator that requires gas flow 100 l/min at the minimum pressure 2.8 bar), check required flow of connecting medical device with pressure and flow characteristics of the valve stated in appendix nr. 1. For assurance of sufficient performance (pressure and flow characteristics of the valve) do not use medical package if pressure indicator enters red zone.**

Flow outlet

- **Before connecting any accessory to the flow outlet make sure that the patient is not connected before using the valve.**

- Ensure that the hose/humidifier is compatible with the flow outlet features.
- Push hose onto the regulator flow outlet/screw-on the humidifier.
- Ensure that the hose/humidifier is well engaged/screwed-on.

6b.3 After use

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Ensure the flow control knob (if any) is on ZERO position.
- Disconnect all downstream equipment from user outlets.
- Refit pressure outlet and flow outlet protection cap. Before refitting the caps, ensure they are clean.

6c. Operations performed by the user / patient

6c.1 Use of user outlet connections

Use of product through the flow outlet (Flow setting)

- Ensure that the flow control knob is on ZERO position
- Ensure that the accessory is connected to the flow outlet.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Set flow control knob on one required of the available flow rates.

- **Always ensure that the flow control knob has been engaged and not placed between two settings otherwise the flow control head will not deliver correct or any flow of medical gas.**

- **The same variants of flow control head can have an “end stop” in between the maximum flow position and the zero position. Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position (during clockwise rotation) or in zero position (anti-clockwise rotation).**

- **Medical gas flow rate must be prescribed by a doctor.**

After completion of the therapy

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Vent gas pressure from downstream equipment.
- Reset flow control knob on ZERO position when it's not venting anymore - until the flow control knob engages
- Disconnect the tube from the flow outlet.

Use of product through pressure outlet

- Ensure that the flow control knob is on ZERO position (if any).
- Ensure that the accessory IS NOT connected to the pressure outlet.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

- **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Connect the accessory to the pressure outlet.

After completion of the therapy

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Vent gas pressure from downstream equipment.
- Disconnect the accessory from the pressure outlet.

6c.2 After use

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.

- Ensure that the flow control knob is on ZERO position (if any).
- Disconnect all downstream equipment from user outlets.
- Refit pressure outlet and flow outlet protection cap. Before refitting the caps, ensure they are clean.

7. Accessory

Bed hanger

- Bed hanger is used with a guard. It is designed for suspension of medical package – cylinder and combination valve with the guard – onto the patient’s bed. The maximum hanger load is with the use of a 5-litre cylinder.

GCE filling adaptors

- **GCE recommends using GCE filling adaptors; Non GCE filling adaptors can damage the non-return valve in the filling port and impact the safety of the product.**

8. Cleaning

Remove dirt with a soft cloth damped in oil free soap water & rinse with clean water. Used cleaning agents shall comply with oxygen cleanliness standards.

Disinfection can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and they are compatible with the product materials (including labels) and gas.

- **Do not use cleaning solutions containing ammonia!**

- **Do not expose to water or any other liquid.**

- **Do not expose to high temperature (such as autoclave).**

9. Service, product life time and maintenance

9.1 Service period and Product Life time

Serial number and manufacture date

A nine-figure serial number stamped on valve body consists of the data as follows:

RR MM XXXXX

RR: year of manufacture

MM: month of manufacture

XXXXX: sequence number of product

For example: Serial number 050300521 indicates the valve manufactured in 2005, in March, with sequence number 521.

Service period

Maximum time in operation of the GCE combination valves before *Overall Maintenance* is 5 years from the manufacturing date.

GCE also recommends that the products removed from gas cylinder undergo *Overall Maintenance*.

It should not be assumed that Pre / Post Filling Checks (refer to Appendix nr 3) and *Overall Maintenance* covers every safety procedure or method required by local regulations, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedures.

Contact GCE for further information on qualification program for customised service period.

LEP – Life Extension Program

LEP – Life Extension Program is qualification program through which GCE is enabled, based on the checks, to extend the use of the valve. The products qualified for “LEP – Life Extension Program“ need not undergo repairs

made after 5 years from the manufacture date of the product and their life time can be extended by one and more years up to the maximum product life time – 10 years (see Appendix nr. 4 “Life Extension Program”).

Product life time and Waste management

Maximum life time of the product is 10 years from the manufacturing date.

At the end of the product’s life time (10 years maximum), the product must be withdrawn from service. The owner of the device shall prevent the reuse of the product and handle the product in compliance with “Directive of European Parliament and Council 2006/12/ES of 5th April 2006 on waste”.

9.2 Maintenance, Repair and Service Operations

Pre/Post filling Checks

Perform functional and leak-tight checks before and after each filling. See Annex 3 for further instructions. Contact GCE if further information is needed.

Overall maintenance

It consists of the preventive activities as follows:

- replacement of critical components,
- cleaning of reused components,
- product retesting.

Overall Maintenance of GCE combination valves must be carried out by GCE. Contact GCE for further information about service available in your area.

Repair operations

Repairs consist of corrective operations (i.e.: involve product known as faulty).

Repair activities of GCE combination valves must be carried out by GCE or authorized repair centres. Contact GCE for further information about service available in your area.

Repair activities can generally be performed with the combination valve still fitted on the cylinder.

Repair activities for all GCE combination valve families, which do not require an authorization for repair, are the replacements of the following damaged components:

- guard,
- labels,
- caps and detachable hose nipples.

Depending on the combination valve type, there are also other replacements/repairs of components, which do not require repair authorisation. Contact GCE for further information.

- **All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

- **Use only genuine GCE components.**

Any product sent back to GCE (or GCE authorised centre) for repair or maintenance shall properly be packaged according GCE recommendations. Contact GCE for further information about available maintenance packaging set. The purpose of the maintenance shall clearly be specified (repair, overall maintenance). For product to be repaired, the fault short description and any reference to a claim nr should be indicated.

10. Glossary

Consult operating instruction

Upper and lower humidity limit

Caution




Keep away from oil and grease

Suitable for Home care use

Keep away from heat and flammable material

 Suitable for Hospital care use	 Product serial number
 Suitable for Emergency care use	 Product part (order) number
 Keep dry	 Inlet parameter
 Upper and lower temperature limit	 Outlet parameter
 Manufacture date	 Manufacturer
 Fragile, handle with care	 Batch number
 The date of service or use The perforated number indicates the year the product has to undergo the overall maintenance activities or has to be removed from service. Refer to the serial number to determine what has to be done	 Life Extension Program
 Atmospheric pressure limitation	 To be used until
 Inlet pressure	 Customer number
	 Outlet pressure

11. Warranty

GCE guarantees pressure regulator one year or in accordance to local regulations, from date of delivery, against faulty design, material & workmanship.

GCE shall not be liable for loss of production, loss of profit or any other consequential damage or indirect loss. In the event of any fault in the goods due to defective design, materials or workmanship, the liability of GCE is limited to replacement of these goods provided that written notification is given to GCE within three months from the date of delivery or deemed delivery or such shorter time as may be specified in the quotation. The goods returned to GCE will be not accepted unless GCE gives previous written consent to their return.

The liability of the GCE Combination valves for the use with medical gases is irrevocably transferred to the owner or operator to the extent that it is modified, serviced or repaired by personnel not employed or authorized by GCE or if the apparatus is used in a manner not conforming to its intended use.

GCE cannot be held responsible for the misuse of the equipment in case of failure to comply with the instructions for use.

APPENDIX

Nr 1 – Technical and performance data

Nr 2 – Quick-coupler features and connecting/disconnecting procedure

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© GCE s.r.o.